
Surgical sequence Implant Ø 3.3 3.75 4.0 5.0

Pilot hole Ø	Ostwill n.	Ostwill n.	Ostwill n.	Implant Ø
1.6 mm	①	②	③	3.3 mm
2.0 mm	②	③	④	3.75 mm
2.0 mm	②	③	⑤	4.0 mm
3.0 mm	⑤	⑥	⑦	5.0 mm

www.metahosp.com

Ref 004172 Ostwill set

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1 autoclavable tray with 7 osteotomes, 7 sterile and disposable depth probes and 4 sterile and disposable torque spanners

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20 sterile blisters containing 7 depth probes and 20 sterile blisters containing 1 white torque spanner and 1 red torque spanner

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Rev 01 07-2005

Ostwill system
Minimally Invasive Osteotomy Technique

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
At last, soft bone will be treated like hard bone

Oswill system, indicated for soft bone treatment, uniquely combines an increase in bone density of the bone-implant contact and an increase in the thickness of the alveolus.

The system is composed of seven re-sterilisable osteotomes, in titanium alloy with a tapered screw design, incremental in diameter and disposable accessories: seven depth probes with markings at 8-10-11.5-13-15 mm and 2 different torque spanners which provide, while screwing, the minimal invasivity on the bone and the atraumaticity for the patient.

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Atraumatic osteotomy

Oswill system does not involve the use of a mallet, which is always traumatic for the patient. The 2 special torque spanners avoid the reaching of high compressive stresses minimizing the risk of vascular damage and bone fractures near the implant site.

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Gradual control of the ridge expansion

The thickness of the ridge increases in a controlled, progressive way in soft bone. The screw tapered osteotome, the specific design of the threads and its incremental effect both compress and expand the implant site.

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Increase in bone density

Bone compression provides a better initial adjustment of the trabecular bone and increased bone-implant contact with an increase in bone density. In vitro tests on sheep trabecular bone show increased bone - implant contact.

(picture on the left side) implant placed after an osteotomy drills.

(picture on the right side) implant placed after *Oswill system*.

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Increase in the primary stability of the implant

There is an increase in the initial stability of the implant during the healing process, resulting in an improvement of the reparative process. During bone compression, the micro-fractured trabeculas are lateralised; this leads to an increase in bone density, followed by intensive, above-normal osteogenetic activity.

Oswill system improves the new minimally invasive techniques for flapless implant placement.

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It can be used with different implant systems

Oswill system can be used with most of the cylindrical and tapered threaded implants.

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Simple and safe surgical technique

Oswill system uses probes, which reduce the screwing depth to 8-10-11.5-13 and 15mm avoiding for the surgeon the risks due to the visual evaluation.

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Maximum control of bone compression

The torque spanners reduce the risk of bone fractures, control the compressive stresses on the bone, possible cause of necrosis, and facilitate the use of the threaded osteotomes in posterior areas. Moreover they allow to check the bone tissue for quality and consistency during the surgery.

The red spanner is indicated when the thickness of the alveolar ridge is more than 6mm or in presence of a particularly hard bone.

Preparazione del foro guida con fresa Ø 2.5 mm
Preparation of the pilot hole with Ø 2.5 mm drill

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